

INSTRUCTIONS FOR USE

IVD

■ 8248 Gastrointestinal Control Panel (Inactivated Swab)

INTENDED USE

The Gastrointestinal Control Panel (Inactivated Swab) is intended for use as non-viable, external, positive and negative control material to evaluate the performance of nucleic acid amplification testing (NAAT) procedures that detect the analytes in Table 1. This product has no qualitative or quantitative assigned value. This control material is nonautomated and not intended to be used for screening, monitoring, or diagnosis. This control is not intended for any specific patient population or specimen.

SUMMARY AND PRINCIPLES

The Gastrointestinal Control Panel (Inactivated Swab) can be used to monitor the extraction, amplification and detection process of molecular testing assays that include the analytes in Table 1. Routine use of quality controls monitors test variation, lot-to-lot test kit performance, operator performance, and aid in identifying random or systemic error.

COMPOSITION

The Gastrointestinal Control Panel (Inactivated Swab) consists of 6 individually packaged lyophilized positive control swabs and 6 individually packaged lyophilized negative control swabs. The analytes in Table 1 have been inactivated using thermal treatments.

The Gastrointestinal Control Panel (Inactivated Swab) is lyophilized in a PCR compatible matrix. The organisms are prepared in a buffered solution with materials of plant and animal origin, preservatives, and stabilizers. The solution is lyophilized into a ready-to-use swab.

Table 1: Contents of the Gastrointestinal Control Panel (Inactivated Swab)

Analytes*
Positive Control
<i>Campylobacter jejuni</i>
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium
<i>Escherichia coli</i> (STEC)
<i>Shigella sonnei</i>
<i>Vibrio cholerae</i>
<i>Vibrio parahaemolyticus</i>
<i>Yersinia enterocolitica</i>
<i>Cryptosporidium parvum</i>
<i>Giardia lamblia</i>
Norovirus GI surrogate
Norovirus GII surrogate
Negative Control
Negative Patient Relevant Matrix**

*All detectable analytes are added at a target concentration of 10^4 - 10^7 copies per swab. These are input concentrations and are not representative of recoverable concentrations or expected values.

***Lactobacillus acidophilus* is added as a component of the negative patient relevant matrix and not intended for positive detection.


HELIX|ELITE™
MOLECULAR STANDARDS


Microbiologics™
OC SETS AND PANELS

WARNINGS AND PRECAUTIONS

- For In Vitro Diagnostic use only.
- For professional use only. To be used by personnel trained in the use of the assay.
- The inactivated lyophilized swabs are single-use only. If hydrated, do not freeze for reuse. If reused, the inactivated lyophilized swabs may not perform as expected including, but not limited to, degradation or over-dilution of material resulting in false negative results, contamination resulting in false positive results or positive detection of organisms not on the swab, and environmental contamination.
- Do not bend or break the swab while it is still within the foil pouch.
- Do not open foil pouch until ready to use.
- Although this product has been inactivated, there is no known test or inactivation method that can assure that it will not transmit infection. This product must be treated as a potential biohazard and be handled using universal laboratory precautions. Wear appropriate personal protective equipment. Do not pipette by mouth. Do not smoke, eat, or drink in areas where specimens are handled. Disinfect any spills and dispose of all materials in accordance with national and local regulations.
- Refer to the Safety Data Sheet (SDS) for more detailed information. The SDS can be located on the Microbiologics website at www.microbiologics.com or by contacting Customer Service at info@microbiologics.com.
- This product does not contain any hazardous substances listed in 1272/2008/EC.
- Report any serious incident that has occurred in relation to the device to Microbiologics and the local regulatory officials in which the user and/or the patient is established.



STORAGE AND EXPIRATION

Store the Gastrointestinal Control Panel (Inactivated Swab) at 2°C-25°C in the original packaging up to the indicated expiration date. After opening the foil pouch, use immediately. Do not use swab more than 1 hour after opening and do not exceed 25°C after pouch has been opened.

The Gastrointestinal Control Panel (Inactivated Swab) should not be used if:

- Stored improperly
- There is evidence of excessive exposure to heat or moisture
- The expiration date has passed
- Packaging or swab is damaged

MATERIALS REQUIRED BUT NOT PROVIDED


- Assay for detection including required instrumentation and consumables

INSTRUCTIONS FOR USE


1. Read package insert, instructions for use or lab protocol for the applicable assay. Some instruments and assays are equipped with special QC settings. In these instances, it may be necessary to use the special setting when using QC sets and panels.
2. Tear open pouch at notch and remove swab from pouch.
3. Insert the swab as instructed by the assay manufacturer. See illustrated instructions on page 5.
 - a. If protocol instructs to break the swab, insert the swab into the device.
 - i. Swab contains an embossed line to facilitate breaking. This breakpoint height may not be appropriate for all assays.
 - ii. If breakpoint is not appropriate for assay, lift the swab up slightly from the bottom of the device and break/snap the swab to the right by pushing the shaft of the swab against the opening of the device. The break point should be appropriate to the size of the test device or vial.
 - iii. Start process over with new swab if not broken/snapped appropriately.
4. Process following assay instructions for use.

LIMITATIONS

- When breaking the control swab, if left too short for the test device or vial, the swab may be jostled and not make sufficient contact with reagents. If swab is left too long, the swab shaft may compromise the integrity of the test device or vial by preventing proper closure. Do NOT cut shaft as contamination may occur.



✘ Swab left too short
Do NOT use.






















✘ Swab left too long. Do NOT use. Do NOT cut shaft to length. Start process over with new swab.

- This product is a DNA control and is not intended as a control for the reverse transcription process.
- There are no known extrinsic factors or interfering substances.
- This product is unassayed control material. It may not be suitable for use with all kits and procedures as not all instruments and assays are compatible with multi-target controls. Customer is responsible for verifying the performance of this product with their chosen instrumentation and assay(s). As a third-party control manufacturer, Microbiologics' provides quality controls that deliver an independent, unbiased assessment of performance with any instrument or method. While not intended to replace control materials provided by the assay/instrument supplier, third-party control materials should be considered.
- Target concentrations of each analyte are specific to Microbiologics' assay method and procedures. These organisms are intact, non-viable, and may be used with any PCR-based test or assay. Microbiologics guarantees each nucleic acid is present and can be amplified but does not guarantee specific analyte concentrations. Each laboratory should establish its own range of acceptable values on their assay system per their internal quality assurance procedure/program. Nucleic acid reactivity, which may vary over time, is dependent on a laboratory's instrumentation, assay method, procedures, calibration, or technician. Microbiologics' molecular controls are not calibrators and should not be used for assay calibration or as an absolute reference material.

MICROBIOLOGICAL STATE

This product was prepared using suitable inactivation methods. While the product has been tested for innocuity, universal laboratory precautions are recommended, and material should be treated as though it was a viable specimen.

KEY OF SYMBOLS

	Authorized Representative in the European Community / European Union		Health hazard
	Batch code (Lot)		In vitro diagnostic medical device
	Biological risks		Manufacturer
	Catalog number		Negative control
	Caution		Positive control
	CE mark		Quantity
	Consult instructions for use or consult electronic instructions for use		Telephone number
	Contains sufficient for <n> tests		Temperature limit
	Do not re-use		Use-by-date
	Do not use if package is damaged and consult instructions for use		

Please refer to product labels for applicable symbols.

NOTICE TO PURCHASERS

The purchase of this product allows the purchaser to use it for Quality Control. No general patents or other license of any kind other than this specific right of use from purchase is granted hereby. No other rights are conveyed expressly, by implication or by estoppel to any other patents. Furthermore, no rights for resale are conferred with the purchase of this product.

The Microbiologics logo is a registered trademark of Microbiologics, Inc.

WEBSITE

Visit our website, www.microbiologics.com, for current technical information and product availability.

ASSISTANCE



Microbiologics, Inc.

200 Cooper Avenue North
St. Cloud, MN 56303 USA
www.microbiologics.com

Technical Support

Tel: +1.320.229.7045
U.S. Toll Free: +1.866.286.6691
Email: techsupport@microbiologics.com

Customer Service


Tel: +1.320.253.7400
U.S. Toll Free: +1.800.599.2847
Email: info@microbiologics.com

Additional copies of this product insert may be obtained at www.microbiologics.com or by emailing info@microbiologics.com

ILLUSTRATED INSTRUCTIONS


Each kit consists of 6 individually packaged lyophilized positive control swabs and 6 individually packaged lyophilized negative control swabs .

1



Read package insert, instructions for use or lab protocol for the applicable assay. Some instruments and assays are equipped with special QC settings. In these instances, it may be necessary to use the special setting when using QC sets and panels.

2




Tear open pouch at notch and remove swab from pouch.

3

Insert the swab as instructed by the assay manufacturer.

If protocol instructs to break the swab, insert the swab into the device. The swab contains an embossed line to facilitate breaking. This breakpoint height may not be appropriate for all assays. If breakpoint is not appropriate for assay, lift the swab up slightly from the bottom of the device and break/snap the swab to the right by pushing the shaft of the swab against the opening of the device. The break point should be appropriate to the size of the test device or vial. Start process over with new swab if not broken/snapped appropriately.



Swab left too short. Do NOT use.

Swab left too long. Do NOT use. Do NOT cut shaft to length. Start process over with new swab.

4

Process following assay instructions for use.

REVISION HISTORY ---

Publication History		
Revision	Date	Description of Change
A	2023-04-27	Initial release.
B	2025-07	Revise the contents of Table 1 and update the wording of the Storage and Expiration section.